



 NovaDose

TRUE
PHARMACEUTICAL

0
1
2

Youth In Your Hands



NovaDose

What Is NovaDose?

pure, unfiltered energy

faster muscle recovery

aesthetic rejuvenation

 NovaDose

NovaDose delivers high-purity, pharmaceutical-grade **NAD+** through a streamlined, user-friendly system that is suitable for both at-home routines and professional clinical use.



True Wellness

CONVENIENCE

Enjoy effortless microdosing wherever you are - it's as easy as click-and-go, with no hassle.

PERSONALIZED WELLNESS

Support energy, focus, and cellular health in a format that fits your lifestyle.



Starts Here

INNOVATION

A modern delivery system offering precise dosing and a smooth, gentle experience.

SAFETY & QUALITY

Manufacturing in GMP-certified facilities to ensure purity consistency, and reliable quality.

Why Choose NovaDose?



PHARMA REGISTERED NAD+

High-quality B-GMP sourced NAD+ for maximum purity and effectiveness.

COST-EFFECTIVE

More affordable and sustainable than IV therapy.

CLINICALLY OPTIMIZED FORMULA

Scientifically dosed for fast absorption & near 100% bioavailability.

Formulation Overview

NAD+

Primary cellular coenzyme supporting mitochondrial energy production and metabolic regulation.

L-CARNITINE


Facilitates fatty acid transport into mitochondria, supporting efficient energy utilization.

TAURINE

Supports cellular balance, membrane stability, and neurological function.

REDUCED GLUTATHIONE

Endogenous antioxidant involved in redox balance and cellular protection.




NovaDose Wellness Benefits


Cellular Longevity



Improved Metabolism




Balanced Mood



Improved Sleep Quality



Enhanced Focus



Brighter Skin



NovaDose - Dosage & Storage

DEVICE FORMAT

- Auto injection pen
- Contains **1000mg NAD+**
- Delivered though 2 cartridges
- Total volume: 6 ml

DOSAGE

- Adjustable dose settings (1-5):
10–100 mg NAD+

Loading Phase

Setting 1 – 10 mg
Every other day

Titration Phase

Setting 2 – 25 mg
Every other day

After titration

Settings 3–5 (50–100 mg)
Adjusted based on individual response

STORAGE

- Store pens and cartridges in the refrigerator
- 2–8°C (35.6–46.4°F)
- Do not freeze
- Keep away from direct sunlight
- Use cartridges before expiration date

Contact Us

Email

info@thetruepharma.com

Phone

+82 10 6399 4897

Korea office

198-26 Hapjeong-dong, Mapo-gu,
Seoul, South Korea

USA office

17654 Newhope St, Ste D, Fountain
Valley, CA, USA

